


LEVITY

by ORA Medical, Inc.

Download Ora Connect for:

- Levity video guides
- Suggested workouts
- Progress tracking

Available on iOS & Android!



Upper/Lower Backplate



- Remove the harness back cushion.



UPPER BACKPLATE

- Securely fasten the three bolts

Tip: Screw all the bolts in partially first, then tighten.



LOWER BACKPLATE

- Securely fasten the three bolts

Tip: Screw all the bolts in partially first, then tighten.



- Use the Velcro strip to attach the cushion to the backplate.

Combined Upper & Lower Backplate



- Align the lower and upper backplate together.

Tip: Place the upper backplate on top.



- Securely fasten the three bolts.

Tip: Screw all the bolts in partially first, then tighten.



- Use the Velcro strip to attach the large cushion to the top backplate.



- Use the Velcro strip to attach the small cushion to the bottom backplate.



Refer to the user manual.

Contact Us

+1 (514) 839-1773
info@oramedical.ca
80 rue Queen, Montréal
QC H3C 2N5
CANADA

Thank you for choosing our product !

The Ora Medical team is excited to accompany you as you take your first steps with Levity!

Share your experience with us on Instagram by tagging @ora_medical.

Chest Harness Installation



- Attach the backplate as per *Upper Backplate Guide*.



- Attach the belt to the harness backplate.




- Attach the chest harness to the buckles on the Upper Backplate.



- Attach the chest harness to the belt.

Chest Harness Adjustment



- Adjust the harness to hip level.
- Place the user in the device, and adjust their pelvic support.



- Attach the belt around the user.

Tip: The belt can be installed on the user prior to installation in the device.



- Pull the chest harness over the user's head and adjust the height.



- Attach the chest harness buckles to the belt.

Tip: User should be upright.



Refer to the user manual.